

# SAMPLES. STORED DRY.

**Provide the optimal environment for your samples with a dry storage cabinet.**

Biomātrica<sup>®</sup> dry storage cabinets employ passive control devices to ensure a stable, dry environment and are useful for laboratories where humidity fluctuations are an issue. The cabinets maintain relative humidity at <40% for the optimal room temperature storage of samples stabilized in DNASTable<sup>®</sup>\*, RNASTable<sup>®</sup> and CloneStable<sup>®</sup> (for bacterial DNA) products.

**Protect your precious samples by utilizing a controlled dry environment.**


- Ensure environmental consistency even in labs with humidity fluctuations.
- Requires no energy input; help your lab save money and “go green.”
- Ideal for use with samples stabilized in Biomātrica technology.

**Select from a variety of configurations to best fit your lab space and work flow.**

**Maintain a constant dry environment at room temperature.**

**Compatible with samples stored in DNASTable<sup>®</sup>, RNASTable<sup>®</sup>, or CloneStable<sup>®</sup>.**

**Configurable to match your lab space and work flow.**



\*Also available as QIAsafe™.

**Maintain a constant dry environment at room temperature.**

- Available in three different sizes- small, medium and large.
- Adjustable internal shelves accommodate a combination of single tubes, multi-well plates and tube-plates.
- Stackable units allow for increased storage capacities.
- Custom configurations are also available.

**Dry storage cabinets are available in the following sizes**

PRODUCT	CATALOG #	SPECIFICATIONS
<p>Small Dry Storage Cabinet</p> 	<p>With (1) rack 95905-058</p>	<p><b>CAPACITY</b> With rack (2): 96w - <b>4,128</b> samples; 0.75ml 2D tube-plates - <b>1,728</b> samples; 384w - <b>16,512</b> samples. Without rack: 96w - <b>6,624</b> samples; 0.75ml 2D tube-plates - <b>2,304</b> samples; 384w - <b>26,496</b> samples</p> <p><b>EXTERNAL DIMENSIONS</b> 8.7" (W) x 11.7" (H) x 10.7" (D)</p>
<p>Medium Dry Storage Cabinet</p> 	<p>With racks 95904-158</p> <p>Without racks 95904-168</p>	<p><b>CAPACITY</b> With racks (6): 96w - <b>14,592</b> samples; 0.75ml 2D tube-plates - <b>5,952</b> samples; 384w - <b>58,368</b> samples. Without racks: 96w - <b>14,592</b> samples; 0.75ml 2D tube-plates - <b>9,216</b> samples; 384w - <b>105,984</b> samples.</p> <p><b>EXTERNAL DIMENSIONS</b> 25" (W) x 11.5" (H) x 10.2" (D)</p>
<p>Large Dry Storage Cabinet</p> 	<p>With racks 95905-158</p> <p>Without racks 95904-178</p>	<p><b>CAPACITY</b> With racks (12): 96w - <b>29,184</b> samples; 0.75ml 2D tube-plates - <b>11,904</b> samples; 384w - <b>116,736</b> samples. Without racks: 96w - <b>29,184</b> samples; 0.75ml 2D tube-plates - <b>18,432</b> samples; 384w - <b>211,968</b> samples.</p> <p><b>EXTERNAL DIMENSIONS</b> 25" (W) x 23" (H) x 10.2" (D)</p>
ACCESSORIES	CATALOG #	DESCRIPTION
Dessicant Pack	17901-080	Sorbent System.
Mini Hygro-Thermometer	16901-090	Digital thermometer to track humidity levels.
Adjustable Cabinet Rack	95904-058	Adjustable rack for dry storage cabinet. Fully adjustable - holds tubes, plates and tube-plates.

To order, please call 866-379-6879 or visit [www.biomatrica.com](http://www.biomatrica.com).